

CAH10 Polyclonal Antibody

Cat No: HR1AP10608

For research use only

Overview

Product Name	CAH10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	CA10 hucep-15 UNQ533/PRO1076
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>carbonic anhydrase 10(CA10) Homo sapiens This gene encodes a protein that belongs to the carbonic anhydrase family of zinc metalloenzymes, which catalyze the reversible hydration of carbon dioxide in various biological processes. The protein encoded by this gene is an acatalytic member of the alpha-carbonic anhydrase subgroup, and it is thought to play a role in the central nervous system, especially in brain development. Multiple transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008].</p>