



CHSTE Polyclonal Antibody

Cat No: HR1AP10606

For research use only

Overview

Product Name	CHSTE Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	CHST14 D4ST1 UNQ1925/PRO4400
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	carbohydrate sulfotransferase 14(CHST14) Homo sapiens This gene encodes a member of the HNK-1 family of sulfotransferases. The encoded protein transfers sulfate to the C-4 hydroxyl of N-acetylgalactosamine residues in dermatan sulfate. Mutations in this gene have been associated with adducted thumb-clubfoot syndrome.[provided by RefSeq, Mar 2010],