

## CHST3 Polyclonal Antibody

Cat No: HR1AP10604

For research use only

### Overview

Product Name	CHST3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	CHST3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	carbohydrate sulfotransferase 3(CHST3) Homo sapiens This gene encodes an enzyme which catalyzes the sulfation of chondroitin, a proteoglycan found in the extracellular matrix and most cells which is involved in cell migration and differentiation. Mutations in this gene are associated with spondylepiphyseal dysplasia and humerospinal dysostosis. [provided by RefSeq, Mar 2009],