

CALL6 Polyclonal Antibody

Cat No: HR1AP10592

For research use only

Overview

Product Name	CALL6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	CALML6 CAGLP CALGP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the calmodulin family. Calglandulin subfamily.,similarity:Contains 4 EF-hand domains.,tissue specificity:Expressed in prostate, thymus, heart, skeleton muscle, bone marrow and ovary.,