

CAN13 Polyclonal Antibody

Cat No: HR1AP10594

For research use only

Overview

Product Name	CAN13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	CAPN13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>calpain 13(CAPN13) Homo sapiens The calpains, calcium-activated neutral proteases, are nonlysosomal, intracellular cysteine proteases. The mammalian calpains include ubiquitous, stomach-specific, and muscle-specific proteins. The ubiquitous enzymes consist of heterodimers with distinct large, catalytic subunits associated with a common small, regulatory subunit. This gene encodes a member of the calpain large subunit family. [provided by RefSeq, Jun 2012].</p>