

CMTA2 Polyclonal Antibody

Cat No: HR1AP10590

For research use only

Overview

Product Name	CMTA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	132kDa
GeneID?Human?	CAMTA2 KIAA0909
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calmodulin binding transcription activator 2(CAMTA2) Homo sapiens The protein encoded by this gene is a member of the calmodulin-binding transcription activator protein family. Members of this family share a common domain structure that consists of a transcription activation domain, a DNA-binding domain, and a calmodulin-binding domain. The encoded protein may be a transcriptional coactivator of genes involved in cardiac growth. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jan 2010],