

CAPS2 Polyclonal Antibody

Cat No: HR1AP10587

For research use only

Overview

Product Name	CAPS2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	142kDa
GeneID?Human?	CADPS2 CAPS2 KIAA1591
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcium dependent secretion activator 2(CADPS2) Homo sapiens This gene encodes a member of the calcium-dependent activator of secretion (CAPS) protein family, which are calcium binding proteins that regulate the exocytosis of synaptic and dense-core vesicles in neurons and neuroendocrine cells. Mutations in this gene may contribute to autism susceptibility. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2009],