

CABP4 Polyclonal Antibody

Cat No: HR1AP10584

For research use only

Overview

Product Name	CABP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	CABP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcium binding protein 4(CABP4) Homo sapiens This gene encodes a member of the CABP family of calcium binding protein characterized by four EF-hand motifs. Mutations in this gene are associated with congenital stationary night blindness type 2B. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2014],