

## CB39L Polyclonal Antibody

Cat No: HR1AP10583

For research use only

### Overview

Product Name	CB39L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CAB39L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	miscellaneous:Found in serum of 50% of patients with acute monocytic leukemia.,similarity:Belongs to the Mo25 family.,