

CHERP Polyclonal Antibody

Cat No: HR1AP10575

For research use only

Overview

Product Name	CHERP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	100kDa
GeneID?Human?	CHERP DAN26 SCAF6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Involved in calcium homeostasis, growth and proliferation.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 CID domain.,similarity:Contains 1 G-patch domain.,similarity:Contains 1 SURP motif repeat.,subcellular location:Distributed throughout the cytoplasm and also localizes to the perinuclear region of both human erythroleukemia (HEL) cells and Jurkat cells. Colocalizes with ITPR1.,tissue specificity:Expressed in brain, placenta, lung, liver, kidney, pancreas, cardiac and skeletal muscle, and in cultured HEL and Dami cells.,</p>