

CABIN Polyclonal Antibody

Cat No: HR1AP10573

For research use only

Overview

Product Name	CABIN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	H,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	244kDa
GeneID?Human?	CABIN1 KIAA0330
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>calcineurin binding protein 1(CABIN1) Homo sapiens Calcineurin plays an important role in the T-cell receptor-mediated signal transduction pathway. The protein encoded by this gene binds specifically to the activated form of calcineurin and inhibits calcineurin-mediated signal transduction. The encoded protein is found in the nucleus and contains a leucine zipper domain as well as several PEST motifs, sequences which confer targeted degradation to those proteins which contain them. Alternative splicing results in multiple transcript variants encoding two different isoforms. [provided by RefSeq, Jan 2011],</p>