

BSSP4 Polyclonal Antibody

Cat No: HR1AP10570

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | BSSP4 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 34kDa |
| GeneID?Human? | PRSS22 BSSP4 PRSS26 SP001LA UNQ302/PRO343 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | protease, serine 22(PRSS22) Homo sapiens This gene encodes a member of the trypsin family of serine proteases. The enzyme is expressed in the airways in a developmentally regulated manner. The gene is part of a cluster of serine protease genes on chromosome 16. [provided by RefSeq, Jul 2008], |