

BRI3 Polyclonal Antibody

Cat No: HR1AP10569

For research use only

Overview

Product Name	BRI3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	BRI3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the BRI3 family.,