

B3GT5 Polyclonal Antibody

Cat No: HR1AP10556

For research use only

Overview

Product Name	B3GT5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	B3GALT5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>beta-1,3-galactosyltransferase 5(B3GALT5) Homo sapiens This gene encodes a member of a family of membrane-bound glycoproteins. The encoded protein may synthesize type 1 Lewis antigens, which are elevated in gastrointestinal and pancreatic cancers. Alternatively spliced transcript variants have been observed for this gene, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2013],</p>