

BNIPL Polyclonal Antibody

Cat No: HR1AP10553

For research use only

Overview

Product Name	BNIPL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	BNIPL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BCL2 interacting protein like(BNIPL) Homo sapiens The protein encoded by this gene interacts with several other proteins, such as BCL2, ARHGAP1, MIF and GFER. It may function as a bridge molecule between BCL2 and ARHGAP1/CDC42 in promoting cell death. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2011],