

IRK15 Polyclonal Antibody

Cat No: HR1AP10549

For research use only

Overview

Product Name	IRK15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	KCNJ15 KCNJ14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>potassium voltage-gated channel subfamily J member 15(KCNJ15) Homo sapiens Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein has a greater tendency to allow potassium to flow into a cell rather than out of a cell. Eight transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2013].</p>