

## RECQ4 Polyclonal Antibody

Cat No: HR1AP10544

For research use only

### Overview

Product Name	RECQ4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	132kDa
GeneID?Human?	RECQL4 RECQ4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>RecQ like helicase 4(RECQL4) Homo sapiens The protein encoded by this gene is a DNA helicase that belongs to the RecQ helicase family. DNA helicases unwind double-stranded DNA into single-stranded DNAs and may modulate chromosome segregation. This gene is predominantly expressed in thymus and testis. Mutations in this gene are associated with Rothmund-Thomson, RAPADILINO and Baller-Gerold syndromes. [provided by RefSeq, Jan 2010],</p>