

ABCF3 Polyclonal Antibody

Cat No: HR1AP10537

For research use only

Overview

Product Name	ABCF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	ABCF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,