

ATPO Polyclonal Antibody

Cat No: HR1AP10527

For research use only

Overview

Product Name	ATPO Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	ATP5O ATPO
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ATP synthase, H⁺ transporting, mitochondrial F1 complex, O subunit(ATP5O) Homo sapiens The protein encoded by this gene is a component of the F-type ATPase found in the mitochondrial matrix. F-type ATPases are composed of a catalytic core and a membrane proton channel. The encoded protein appears to be part of the connector linking these two components and may be involved in transmission of conformational changes or proton conductance. [provided by RefSeq, Jul 2008].</p>