



ARSJ Polyclonal Antibody

Cat No: HR1AP10518

For research use only

Overview

Product Name	ARSJ Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	ARSJ UNQ372/PRO708
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	arylsulfatase family member J(ARSJ) Homo sapiens Sulfatases (EC 3.1.5.6), such as ARSJ, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules (Sardiello et al., 2005 [PubMed 16174644]).[supplied by OMIM, Mar 2008],