



## AAA1 Polyclonal Antibody

Cat No: HR1AP10519

For research use only

### Overview

Product Name	AAA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	SLC7A10 ASC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 7 member 10(SLC7A10) Homo sapiens SLC7A10, in association with 4F2HC (SLC3A2; MIM 158070), mediates high-affinity transport of D-serine and several other neutral amino acids (Nakauchi et al., 2000 [PubMed 10863037]).[supplied by OMIM, Mar 2008],