

AQP6 Polyclonal Antibody

Cat No: HR1AP10512

For research use only

Overview

Product Name	AQP6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	AQP6 AQP2L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>aquaporin 6(AQP6) Homo sapiens The protein encoded by this gene is an aquaporin protein, which functions as a water channel in cells. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). This protein is specific for the kidney. This gene and related family members AQP0, AQP2, and AQP5 reside in a cluster on chromosome 12q13. [provided by RefSeq, Jul 2008],</p>