

# AQP10 Polyclonal Antibody

Cat No: HR1AP10511

For research use only

## Overview

Product Name	AQP10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	AQP10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>aquaporin 10(AQP10) Homo sapiens This gene encodes a member of the aquaglyceroporin family of integral membrane proteins. Members of this family function as water-permeable channels in the epithelia of organs that absorb and excrete water. This protein was shown to function as a water-selective channel, and could also permeate neutral solutes such as glycerol and urea. [provided by RefSeq, Jul 2008],</p>