

## B3A4 Polyclonal Antibody

Cat No: HR1AP10506

For research use only

### Overview

Product Name	B3A4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	108kDa
GeneID?Human?	SLC4A9 AE4 SBC5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 4 member 9(SLC4A9) Homo sapiens The protein encoded by this gene is a membrane protein involved in anion exchange. Expression of this gene is mostly limited to the kidney. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012],