

# ANGP4 Polyclonal Antibody

Cat No: HR1AP10502

For research use only

## Overview

Product Name	ANGP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	ANGPT4 ANG3 ANG4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>angiopoietin 4(ANGPT4) Homo sapiens Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The mechanism by which they contribute to angiogenesis is thought to involve regulation of endothelial cell interactions with supporting perivascular cells. The protein encoded by this gene functions as an agonist and is an angiopoietin. [provided by RefSeq, Jul 2008],</p>