

B3A3 Polyclonal Antibody

Cat No: HR1AP10505

For research use only

Overview

Product Name	B3A3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	135kDa
GeneID?Human?	SLC4A3 AE3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 4 member 3(SLC4A3) Homo sapiens The protein encoded by this gene is a plasma membrane anion exchange protein. The encoded protein has been found in brain, heart, kidney, small intestine, and lung. [provided by RefSeq, May 2016],