

# ANGL5 Polyclonal Antibody

Cat No: HR1AP10503

For research use only

## Overview

Product Name	ANGL5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	ANGPTL5 UNQ5795/PRO19600
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 fibrinogen C-terminal domain.,tissue specificity:Mainly expressed in adult heart.,