

PDS5B Polyclonal Antibody

Cat No: HR1AP10499

For research use only

Overview

Product Name	PDS5B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	159kDa
GeneID?Human?	PDS5B APRIN AS3 KIAA0979
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PDS5 cohesin associated factor B(PDS5B) Homo sapiens This gene encodes a protein that interacts with the conserved protein complex termed cohesin. The cohesin complex holds together sister chromatids and facilitates accurate chromosome segregation during mitosis and meiosis. This protein is also a negative regulator of cell proliferation and may be a tumor-suppressor gene. [provided by RefSeq, Jul 2015],