

CTNL1 Polyclonal Antibody

Cat No: HR1AP10492

For research use only

Overview

Product Name	CTNL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	CTNNAL1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:May modulate the Rho pathway signaling by providing a scaffold for the Lbc Rho guanine nucleotide exchange factor (ARHGEF1).,induction:Down-regulated in cancer pancreatic cells undergoing differentiation and apoptosis.,similarity:Belongs to the vinculin/alpha-catenin family.,subunit:Interacts with ARHGEF1.,tissue specificity:Widely expressed. Expressed at lower level in neural tissues and at the highest level in the adrenal gland.,</p>