

A1ATR Polyclonal Antibody

Cat No: HR1AP10490

For research use only

Overview

Product Name	A1ATR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	SERPINA2 ARG5 ATR PIL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	serpin family A member 2 (gene/pseudogene)(SERPINA2) Homo sapiens This gene encodes a member of the serpin family of proteins, a group of proteins that inhibit serine proteases. The encoded intracellular glycoprotein is localized at the endoplasmic reticulum. This gene is a polymorphic pseudogene, with the non-functional allele being predominant in some populations. Some individuals, as represented by the reference genome allele, contain a 2kb coding region deletion and a start code mutation. [provided by RefSeq, Feb 2014].