

# A1AG1 Polyclonal Antibody

Cat No: HR1AP10488

For research use only

## Overview

Product Name	A1AG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	ORM1 AGP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	orosomucoid 1(ORM1) Homo sapiens This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression. [provided by RefSeq, Jul 2008],