

AL1A3 Polyclonal Antibody

Cat No: HR1AP10484

For research use only

Overview

Product Name	AL1A3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	ALDH1A3 ALDH6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	aldehyde dehydrogenase 1 family member A3(ALDH1A3) Homo sapiens This gene encodes an aldehyde dehydrogenase enzyme that uses retinal as a substrate. Mutations in this gene have been associated with microphthalmia, isolated 8, and expression changes have also been detected in tumor cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014],