

AL3B2 Polyclonal Antibody

Cat No: HR1AP10485

For research use only

Overview

Product Name	AL3B2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	ALDH3B2 ALDH8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>aldehyde dehydrogenase 3 family member B2(ALDH3B2) Homo sapiens This gene encodes a member of the aldehyde dehydrogenase family, a group of isozymes that may play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. The gene of this particular family member is over 10 kb in length. The expression of these transcripts is restricted to the salivary gland among the human tissues examined. Alternate transcriptional splice variants have been characterized. [provided by RefSeq, Jul 2008].</p>