

RP25L Polyclonal Antibody

Cat No: HR1AP10482

For research use only

Overview

Product Name	RP25L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	RPP25L C9orf23
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ribonuclease P/MRP subunit p25 like(RPP25L) Homo sapiens This gene encodes a protein that appears to belong to a family of evolutionarily related proteins (DUF78), that may share one or more domains in common. Members of this family are small archaeobacterial proteins with no known function. Alternative splicing has been observed at this locus and two variants, both encoding the same protein, have been identified. [provided by RefSeq, Jul 2008],