

AT2L1 Polyclonal Antibody

Cat No: HR1AP10481

For research use only

Overview

Product Name	AT2L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	AGXT2L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cofactor:Pyridoxal phosphate .,similarity:Belongs to the class-III pyridoxal-phosphate-dependent aminotransferase family.,subunit:Homotetramer.,