

AFF4 Polyclonal Antibody

Cat No: HR1AP10478

For research use only

Overview

Product Name	AFF4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	127kDa
GeneID?Human?	AFF4 AF5Q31 MCEF HSPC092
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>AF4/FMR2 family member 4(AFF4) Homo sapiens The protein encoded by this gene belongs to the AF4 family of transcription factors involved in leukemia. It is a component of the positive transcription elongation factor b (P-TEFb) complex. A chromosomal translocation involving this gene and MLL gene on chromosome 11 is found in infant acute lymphoblastic leukemia with ins(5;11)(q31;q31q23). [provided by RefSeq, Oct 2011],</p>