

ATS20 Polyclonal Antibody

Cat No: HR1AP10472

For research use only

Overview

Product Name	ATS20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	210kDa
GeneID?Human?	ADAMTS20
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ADAM metalloproteinase with thrombospondin type 1 motif 20(ADAMTS20) Homo sapiens The protein encoded by this gene is a member of the ADAMTS family of zinc-dependent proteases. The encoded protein has a signal peptide that is cleaved to release the mature peptide, which is secreted and found in the extracellular matrix. This protein may be involved in tissue remodeling. [provided by RefSeq, Sep 2011],