

ATS13 Polyclonal Antibody

Cat No: HR1AP10470

For research use only

Overview

Product Name	ATS13 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	156kDa
GeneID?Human?	ADAMTS13 C9orf8 UNQ6102/PRO20085
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ADAM metalloproteinase with thrombospondin type 1 motif 13(ADAMTS13) Homo sapiens This gene encodes a member of a family of proteins containing several distinct regions, including a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. The enzyme encoded by this gene specifically cleaves von Willebrand Factor (vWF). Defects in this gene are associated with thrombotic thrombocytopenic purpura. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],</p>