

ADA30 Polyclonal Antibody

Cat No: HR1AP10466

For research use only

Overview

Product Name	ADA30 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	ADAM30 UNQ2509/PRO5997
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ADAM metallopeptidase domain 30(ADAM30) Homo sapiens This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. This gene is testis-specific and contains a polymorphic region, resulting in isoforms with varying numbers of C-terminal repeats. [provided by RefSeq, Jul 2008],</p>