

ADA23 Polyclonal Antibody

Cat No: HR1AP10463

For research use only

Overview

Product Name	ADA23 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	91kDa
GeneID?Human?	ADAM23 MDC3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ADAM metallopeptidase domain 23(ADAM23) Homo sapiens This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. It is reported that inactivation of this gene is associated with tumorigenesis in human cancers. [provided by RefSeq, May 2013],</p>