

# LYPA2 Polyclonal Antibody

Cat No: HR1AP10461

For research use only

## Overview

Product Name	LYPA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	LYPLA2 APT2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lysophospholipase II(LYPLA2) Homo sapiens Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. There are alternatively spliced transcript variants described for this gene but the full length nature is not known yet. [provided by RefSeq, Jul 2008],