

AOAH Polyclonal Antibody

Cat No: HR1AP10460

For research use only

Overview

Product Name	AOAH Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	AOAH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>acyloxyacyl hydrolase(AOAH) Homo sapiens This locus encodes both the light and heavy subunits of acyloxyacyl hydrolase. The encoded enzyme catalyzes the hydrolysis of acyloxylacyl-linked fatty acyl chains from bacterial lipopolysaccharides, effectively detoxifying these molecules. The encoded protein may play a role in modulating host inflammatory response to gram-negative bacteria. Alternatively spliced transcript variants have been described.[provided by RefSeq, Apr 2010],</p>