

ACBD4 Polyclonal Antibody

Cat No: HR1AP10455

For research use only

Overview

Product Name	ACBD4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	ACBD4 HMFT0700
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	acyl-CoA binding domain containing 4(ACBD4) Homo sapiens This gene encodes a member of the acyl-coenzyme A binding domain containing protein family. All family members contain the conserved acyl-Coenzyme A binding domain, which binds acyl-CoA thiol esters. They are thought to play roles in acyl-CoA dependent lipid metabolism. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008],