

ABHDD Polyclonal Antibody

Cat No: HR1AP10446

For research use only

Overview

Product Name	ABHDD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	ABHD13 C13orf6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,