

ABHD3 Polyclonal Antibody

Cat No: HR1AP10447

For research use only

Overview

Product Name	ABHD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	ABHD3 LABH3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>abhydrolase domain containing 3(ABHD3) Homo sapiens This gene encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a very wide range of enzymes. The function of this protein has not been determined. [provided by RefSeq, Jul 2008],</p>