

ABHDA Polyclonal Antibody

Cat No: HR1AP10445

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | ABHDA Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 33kDa |
| GeneID?Human? | ABHD10 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | <p>abhydrolase domain containing 10(ABHD10) Homo sapiens This gene encodes a mitochondrially-localized enzyme that acts in liver cells as a hydrolase. The encoded protein removes glucuronide from mycophenolic acid acyl-glucuronide. There is a pseudogene for this gene on chromosome 6. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013],</p> |