

ABHD1 Polyclonal Antibody

Cat No: HR1AP10444

For research use only

Overview

Product Name	ABHD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	ABHD1 LABH1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>abhydrolase domain containing 1(ABHD1) Homo sapiens This gene is a member of the AB hydrolase superfamily and encodes a protein with an alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions. [provided by RefSeq, Jul 2008],</p>