

# NIP7 Polyclonal Antibody

Cat No: HR1AP10441

For research use only

## Overview

Product Name	NIP7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	NIP7 CGI-37 HSPC031 HSPC180 OK/SW-cl.76 OK/SW-cl.78
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Required for proper 27S pre-rRNA processing and 60S ribosome subunit assembly.,similarity:Belongs to the NIP7 family.,similarity:Contains 1 PUA domain.,subunit:Monomer. Interacts with pre-ribosome complex. May bind to RNA (By similarity). Interacts with NOL8. May interact with SBDS.,