

RS4Y2 Polyclonal Antibody

Cat No: HR1AP10426

For research use only

Overview

Product Name	RS4Y2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	RPS4Y2 RPS4Y2P
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ribosomal protein S4, Y-linked 2(RPS4Y2) Homo sapiens The protein encoded by this gene is a ribosomal protein that is highly similar to RPS4Y1. This gene is located in the male-specific region of the Y chromosome. [provided by RefSeq, Aug 2012],