

# RS26L Polyclonal Antibody

Cat No: HR1AP10424

For research use only

## Overview

Product Name	RS26L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	RPS26P11 RPS26L1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:Could be the product of a pseudogene.,similarity:Belongs to the ribosomal protein S26e family.,